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## Bibliography:

<sup>1)</sup>Health Risk Assessment of VOCs from Surgical Smoke

Shaohua She, Gang Lu, Wah Yang, Mianwei Hong, and Lingfei Zhu; Creative Commons CC; 7/2017

<sup>2)</sup>The Examination of Problems Experienced by Nurses and Doctors Associated with Exposure to Surgical Smoke and the Necessary Precautions;

Journal of Clinical Nursing; 7/2016.

<sup>3)</sup>Tackle Toxic Surgery Smoke Risk for Doctors, Nurses, Patients

John Weekes; Stuff.co.nz; 5/2016.

<sup>4)</sup>Study of Health Care Workers Shows Prevalence of Surgical Smoke Exposure

Andrea Steege; Safety and Health Magazine; 11/2015.

<sup>5)</sup>Surgical Smoke and the Dermatologist

Jesse Lewin, Jeremy Brauer, Ariel Ostad; American Academy of Dermatology, Inc., 2010, doi:10.1016/j.jaad.11.017; 1/2010.

<sup>6)</sup>Surgical Smoke and Infection Control

E. Alp, D.Bijl; R.P. Bleichrodt; B. Hansson; A. Voss; Journal of Hospital Infection; 1/2006.

<sup>7)</sup>BViral Disease Transmitted by Laser-Generated Plume (Aerosol)

Jerome M. Garden, MD; M.Kerry O'Banion, MD; PhD, Abnoeal D. Bakus, PhD; Carl Olson, DVM, PhD; Archives of Dermatology, Vol 138; 10/2002.

<sup>8)</sup>HPV Positive Tonsillar Cancer in Two Laser Surgeons: Case Reports

Margo Rioux, Andrea Garland, Duncan Webster and Edward Reardon; Journal of Otolaryngology, Head and Neck Surgery, 42:54; 1/2013.

<sup>9)</sup>OMaximizing Surgeon Safety During Excimer Laser Photorefractive Procedures

Jay S. Pepose, MD, PhD; American Journal of Ophthalmology, Vol 130, No 6; 12/2002.

<sup>10)</sup>Pilot Study of Directional Airflow and Containment of Airborne Particles in the Size of Mycobacterium Tuberculosis in an Operating Room

Russell N. Olmsted, MPH, CIC; American Journal of Infection Control, Vol 36, No 4; 5/2008.

<sup>11)</sup>Surgical Smoke and Ultrafine Particles

Irene Bröske-Hohlfeld; Gerhard Preissler; Karl-Walter Jauch; Mike Pitz; Dennis Nowak; Annette Peters; H-Erich Wichmann; Journal of Occupational Medicine and Toxicology; 1/2008.

<sup>12)</sup>Respirable Particles in the Excimer Laser Plume

Michael J. Taravella, MD; EyeWorld, Vol 6, No 6; 6/2001.

# Regulations, Recommended Practices & Standards

## Occupational Safety & Health Administration (OSHA)

OSHA is the only U.S. regulatory body to date that has legal authority in the United States granted by Congress. They estimate that 500,000 healthcare workers are exposed to surgical smoke and bio-aerosols each year. On a number of occasions OSHA has reiterated that the management of surgical plume is a healthcare worker safety issue.

### **Association of periOperative Registered Nurses (AORN)**

AORN is a professional association based in Denver, Colorado that represents the interests of more than 160,000 perioperative nurses.

The following are contained in Recommended Practices:

#### **2017 AORN Guideline for Surgical Smoke Safety**

##### **Recommendation I**

"The health care organization should provide a surgical smoke free environment."

##### **Recommendation II**

The perioperative team should evacuate all surgical smoke. "The collective evidence, standards, and guidelines from NIOSH, the Healthcare Infection Control Practices Advisory Committee, and professional organizations indicates that evacuating surgical smoke protects patients and health care workers from the hazards of surgical smoke."

##### **Recommendation III**

"Perioperative team members should receive initial and ongoing education and competency verification on surgical smoke safety."

##### **Recommendation IV**

"Policies and procedures for surgical smoke safety should be developed, reviewed periodically, revised as necessary, and readily available in the practice setting in which they are used."

##### **Recommendation V**

"Perioperative personnel should participate in a variety of quality assurance and performance activities that are consistent with the health care organization's plan to improve understanding and compliance with the principles and processes of surgical smoke safety."

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**HAZARD**

During surgical procedures using a laser or electrothermal unit, the thermal destruction of tissue creates a smoke byproduct. Research studies have confirmed that this smoke plume can contain toxic gases and vapors such as benzene, hydrogen cyanide, and formaldehyde, hexamethyldisiloxane, and fine cellular material including blood fragments, and spores. At high concentrations, the smoke causes visible and upper respiratory tract irritation to health care personnel, and creates visual problems for the surgeon. The smoke has significant odors and has been shown to have mutagenic potential.

**CONTROLS**

NIOSH research has shown airborne contaminants generated by these surgical devices can be effectively controlled. Two methods of control are recommended:

- **VENTILATION**

Recommended ventilation techniques include a combination of general room and local exhaust ventilation (LEV). General room ventilation is not by itself sufficient to capture contaminants generated at the source. The two major LEV approaches used to reduce surgical smoke levels for health care personnel are portable smoke evacuators and room suction systems.

Smoke evacuators contain a suction unit (vacuum pump), filter, hose, and an inlet nozzle. The smoke evacuator should have high efficiency in airborne particle reduction and should be used in accordance with the manufacturer's recommendations to achieve maximum efficiency. A capture velocity of about 100 to 150 feet per minute at the inlet nozzle is generally recommended. It is also important to choose a filter that is effective in collecting the contaminants. A High Efficiency Particulate Air (HEPA) filter or equivalent is recommended for trapping particulates. Various filtering and cleaning processes also exist which remove or inactivate airborne gases and vapors. The various filters and absorbers used in smoke evacuators require monitoring and replacement on a regular basis and are considered a possible biohazard requiring proper disposal.

Room suction systems can pull at a much lower rate and were designed primarily to capture liquids rather than particulates or gases. If these systems are used to capture generated smoke, users must install appropriate filters in the line, ensure that the line is cleared, and that filters are disposed properly. Generally speaking, the use of smoke evacuators are more effective than room suction systems to control the generated smoke from nonendoscopic laser/electric surgical procedures.

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## Australian College of Operating Room Nurses (ACORN)

### Standard S20

- Personnel shall utilize appropriate equipment and procedures to prevent exposure to surgical plume
- Exposure to surgical plume shall be minimized during the surgical procedure
- Surgical smoke capture devices shall be available for use during procedures in which surgical smoke is generated (ACORN 2006)

## Canadian Standards Association (CSA)

CSA Z305.13-13 Plume Scavenging In Surgical, Diagnostic, Therapeutic, and Aesthetic Settings

- Facility policies and procedures shall be written in accordance with (IAW) this Standard.
- Plume shall be evacuated IAW this standard.
- If a facility employs techniques that create plume, they shall have policies that address the potential hazards.

## Danish Working Environment Authority

AT-Instructions 4/2007 and 11/2008

- It is mandatory to implement a measurable setup for local evacuation of harmful substances, such as surgical smoke.
- Such a setup must be equipped with a monitoring feature to indicate if the evacuation system's suction is inadequate.
- Surgical smoke should be removed with local evacuation and as close to the source as possible.
- The filtered air must lead out into the open (read: outside the OR).

## Medicines and Healthcare Products Regulatory Agency (MHRA)

MHRA DB2008(03) April 2008

- Recommends that smoke evacuation systems are to be used during laser surgery. In addition, it is specified that masks and operating room laminar flow systems are not suitable for protection from surgical smoke.

## Association for Perioperative Practice (AfPP)

Standard 2.6 Lasers - Standards and Recommendations for Safe Perioperative Practice

- States that 'Dedicated smoke evacuation machines must be used to remove the smoke...' (AfPP 2007).AT-Instructions 4/2007 and 11/2008
- It is mandatory to implement a measurable setup for local evacuation of harmful substances, such as surgical smoke.

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